

# BeamLine Operations and Safety Awareness (BLOSA) Checklist

# Beamline X28B

All users must be instructed in operating the beamline safely. Leave checkbox blank if not applicable. Training valid 2 years. Visitors use Visitor/Escort forms.

### ACCESS

- Diagnostics Beamline B This is a "Diagnostic & Instrumentation" beamline. Only Machine Operators have access; no general users.

### EMERGENCIES

- Exits R Locate routes to nearest exits
- Fire Extinguishers R Locate fire extinguisher
- Alarm Pulls R Locate fire alarm pull
- Eye Wash/Shower R Locate eye wash/shower [For labs without eye wash stations, prop open door if using corrosives]
- Spill Station R Locate spill control station
- Green/Yellow Boards R Locate/discuss Green Board (beamline emergency contacts, phones and info) and Yellow Board (safety info)

### CONTACTS

- OPCO R Refer to instructions posted on the beamline phone for Operations Coordinator (OPCO) Assistance

### BEAMLINE OPERATIONS

- Enabling Beamline R Review procedures to enable beamline
- Red & Yellow Tags R Provide information about any beamline equipment or systems that are yellow or red tagged
- Beamline Config R Configuration changes to beamline are to be completed by Beamline Staff only
- Power Failure A Identify circuit breaker location; contact OPCOs to reset circuit breakers
- Water Valves A User may adjust water shut-off valves (review location and procedure)
- Water Valves B Water valves are adjusted remotely.
- Air Valves A User may adjust air shut-off valves but requires approval by Beamline Staff (review location and procedure)
- Pink Cards R Beamline may be unattended up to 24 hours unless SAF states "no unattended operations" (review procedure)

### END STATION / EQUIPMENT

- Computer Ops R Reminder to review computer operations, control software, data acquisition software
- Optics A Review optical system set up and procedures
- Optics B Optics (camera) are operated remotely.

### DOCUMENTATION

- Beamline Manuals R Review location of manuals and beamline documentation (must be readily available)
- Reminders R Identify location of Experimental Reminders list posted at the beamline

### EXPERIMENTAL HAZARDS

- Radiation Hazards R Identify radiation locations; inform user to move away from area and call control room if chipmunks sound off
- Radiation UV Ring A Move away from UV ring, mezzanine, platform, transfer lines, and front end during injection
- Electrical Work A No work on exposed electrical components >50V
- Vent System A Review use and procedures for ventilation system and for ventilation alarms, if present
- Hutch Door B Use caution when opening/closing door (tight space, heavy door, pinch hazard)
- Hazard Analysis R Identify location of beamline hazard analysis form

I understand the instructions given to me on beamline operations and safety awareness.

Designated BLOSA trainers for this beamline:

PRINT User Name	Guest #	User Signature	Date	Trainer's Signature	√ UAdm

- Steve Kramer
- Xi Yang