

X17A SAFETY CHECKLIST			BROOKHAVEN NATIONAL LABORATORY NATIONAL SYNCHROTRON LIGHT SOURCE
Rev: A	Effective: 10/24/2010	Page 1 of 1	Number: LS-SCL-0009
Prepared by: Philip Marino		Reviewed by: P. Marino	Approved by: B. Chmiel

Original signatures on file.

The only official copy of this file is the one on-line in the NSLS Quality Assurance website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS QA website.

- 1. Valid Padlock Index and log, no locks open
- 2. Cable penetration shield upstream of double plug doors in place as per photo.
- 3. Cement bricks in place along entire bottom length of the two X17 front end sawtooth plug doors (see photo)
- 4. Labeled Leaded Glass on hutch door.
- 5. Beam stop

OPCO / Beamline Staff: _____ Date: _____

USER SYSTEM CHECKS

- A. Check the status on the PLC for water.
- B. Check the status on the PLC for vacuum.

Beamline Rep.: _____ Date: _____