



# Combined In Situ X-ray Absorption and Diffuse Reflectance Infrared Spectroscopy: A Novel Tool for Catalytic Investigations

N. Marinkovic<sup>1</sup>, Q. Wang<sup>3</sup>, L. Barrio<sup>2</sup>, S.N. Ehrlich<sup>1</sup>, C. Cooper<sup>3</sup>, A.I. Frenkel<sup>3</sup>

<sup>1</sup>University of Delaware; <sup>2</sup>Brookhaven National Laboratory; <sup>3</sup>Yeshiva University

- Investigating the mechanism of catalysis requires complementary spectroscopic techniques. Vibrational spectroscopies are often a method of choice to probe metal-adsorbate interactions. X-ray spectroscopy (XAS) is used to study bulk structure of the catalyst. While diffuse reflectance Fourier transform infrared spectroscopy (DRIFTS) is commonly used on a bench top, XAS requires synchrotron light. To merge the techniques for in-situ catalysis studies, an XAS /DRIFTS instrument was developed at NSLS beamline X18A.
- A rapid-scan FTIR spectrometer is outfitted with a range of detectors and beamsplitters for both mid- and far-infrared measurements, and is further equipped with a modified praying mantis DRIFTS attachment and a custom-made in situ reaction cell. The reaction cell can be heated to 800 °C, allows passage of gases through the catalyst, and is equipped with glassy carbon windows for x-ray light passage. Thus, both XAS and DRIFTS techniques can be done simultaneously with controlled conditions. Along with fast-moving monochromator for quick-EXAFS and mass-spectrometric residual gas analysis, this new tool is a powerful method for testing catalytic reactions in real time.



Experimental setup (left) and the reaction cell (insert in top right). CO adsorption during the reduction of 5% Pt/ $\gamma$ -Al<sub>2</sub>O<sub>3</sub> at elevated temperatures: as the Pt is reduced, which is monitored by XANES at the Pt L<sub>3</sub> edge (middle right), CO is observed to be adsorbed on Pt/ $\gamma$ -Al<sub>2</sub>O<sub>3</sub> by DRIFTS (bottom right)